

Programing the spinners

- **Using the interface panel**
- **View / edit recipes**
- **Operating a program**

- Before programing the spinner, check the correct RPM for the PR you use and its required height. **Use the PR manufacturer's web sites.**

Please note:

- **DO NOT do initialization!**
- **DO NOT change programs other users programs!**

Using the interface panel



Interface panel

[ESC] : move back; previous menu; quit error message

[E/R] : move forward to next page / menu; select and edit process parameter

[Q/I] : confirm error message

↑ / ↓ : move forward / backward

Use the arrows to select the required option and press [E/R].

To move to the next screen select “next” (at the right bottom corner of the screen) and press [E/R].

When a parameter is blinking, press [E/R] first.

View / edit recipes

- Using the arrow, select “Recipe editor” and press [E/R].
- Choose the required program (DO NOT change programs which are not yours!).
- The interface panel will present the program’s parameters.
- After each parameter there’s an option to save the changes or to exit without saving. Choose the required option.
- Refer to “*Main instructions for Spinner*” doc. for proper program parameters and other important rules regarding spinner handling, e.g., PR pouring, etc.

Operating a program

Spinner options

	<p>“F”: when no time was set for a step in a program, “F” will re-activate the program and will move it to the next step.</p> <p>“VAC” : Switch chuck vacuum on / off.</p> <p>“I” : start a program.</p> <p>“O” :Stop a program:</p>
---	--

- When your program is ready, place the wafer on the chuck (refer to “*Placing chucks*” doc.)

Then:

1. Press "VAC".

Interface panel

2. Go to "Automatic" and press [E/R].
3. Select the required program and press [E/R].
4. Go to next and press [E/R]

5. Pour the PR
6. Press "I".

At the end of the program the vacuum is released automatically.

For any inquire call Avraham 052 2291 270